


XZ™ Nanoliter Dialysis Crystallography Plate

The XZ™ plate is a dialysis-based microfluidic protein crystallization device that uses patent-pending microfluidic technology to reduce the consumption of precious protein samples by more than 3,000 times the volumes of the existing dialysis protein crystallization products.

Product Highlights

- Protein consumption as low as 15nL per dialysis chamber
- Fast set-up in as little as 4 minutes for up to 96 dialysis chambers
- Low capital equipment costs to begin running experiments
- High-throughput compatibility with dispensing and imaging systems
- Clear path to structure due to directly scalable designs, easy crystal retrievability and in-situ diffraction capability



Exploded View of the Microfluidic Dialysis Plate

Proven Crystallization Method

Dialysis protein crystallization is a well-established method for protein crystallization. This method utilizes the diffusion and equilibration of small precipitant molecules through a semi-permeable membrane as a means of slowly approaching the concentration at which the protein solute crystallizes.

Precise Nanoliter Metering of Protein Samples

The XZ™ plate patent-pending microfluidic technology enables precise loading of protein samples. Dialysis protein crystallization experiments having protein consumption as low as 15nL per chamber make the XZ™ plate an attractive option for screening experiments.

Directly Scalable Design

The larger volume XZ™ optimization and growth plates (up to 1μL per chamber) have the same diffusion lengths (400μm) to that of the screening plates, therefore, these designs allow for straight forward scale-up to run optimization experiments as well as for translating screening conditions to grow harvestable crystals for use in x-ray diffraction.

Easy Crystal Extraction

Unlike other enclosed microfluidic systems where crystals are difficult to extract, crystals grown on XZ™ plates can be easily extracted because the dialysis membrane in each reagent well can be conveniently removed to allow easy loop access for crystal manipulation.

In-Situ Diffraction Capability

The XZ™ plates can be used as is for meaningful in-situ diffraction data collection using a SBS plate goniometer (BL8.3.1, ALS), thus avoiding time consuming crystal manipulation steps and damages associated with crystal handling.

Availability

1x96 XZ™ Screening, Optimization and Growth Plates are now available. Other configurations are available upon request.


Contact

To start running your crystallization experiments using the XZ™ plates please contact sales@gnbiosystems.com.
+1-510-333-3488. www.gnbiosystems.com


15nL dialysis chamber – Screening Plate
90x




350nL dialysis chamber – Optimization Plate
20x



Lysozyme crystals in screening and optimization dialysis plates



extracted crystal



1.5Å resolution in-situ diffraction data